

FIG. 1




FIG. 2




FIG.3




FIG. 4




FIG. 5




FIG. 6




FIG. 7

STORAGE INTERCHANGE FORMAT




FIG. 8





FIG. 9

[illegible]




FIG. 11




FIG. 12

120

The screenshot shows a DICOM Event Viewer window titled "Event". The window contains a tree view on the left with the following structure:

- Event
 Asir
 Patient
 Birthdate: 1930 9 9
 ID: P0104
 Name: Christiansson, B.
 Affected SOF Instance UID: 1.2.528.1.1001.1.19961112013.200.4.19961117.162655921
 Study
 ID: S023
 Dicom UID: 1.2.528.1.1001.1.960113006.200.2.19960305.183737397
 Sex: M
 S
 R
 A
 T
 D
 S
 Body Part: head
 Manufacturing Equipment
 Modality: US
 Image
 Private File:
 Pixel Rep: 0
 High Bit: 7
 Bits Stored: 8

An "Edit Tag Value" dialog box is open, showing the "Tag Value" field with the text "1930 9 9". The dialog box has "OK" and "Cancel" buttons. A "122" label points to this dialog box.

FIG. 13

124

Patient Info		Additional Patient Info	
ID:	P0104	Other Patient IDs:	P0103
General Patient Info		Other Patient Names:	John Doe
Dicom Family Name Complex:	Christianson, Jason B. Dr. M.D.	Medical Alerts:	None
Prefix:	Dr.	Contrast Allergies:	Not available
First Name:	Jason	Additional History:	
Middle Name:	B.	Yes	
Last Name:	CHRISTIANSSON	Comments:	
Suffix:	MD	Test	
Birth Name:	Smith		
Birthday:	9/9/1930		
Birth Time:			
Mother's Birth Name:	Johnson		
		OK	Cancel

FIG. 14

126

The screenshot shows a dialog box titled "Image Info" with a close button (X) in the top right corner. The dialog contains several fields for image metadata, organized into two columns. The left column includes fields for SOP Class UID, SOP Instance UID, Number, Overlay Number, Lut Number, and Curve Number. The right column includes fields for Key Image, Patient Orientation, Time, Date, and Series UID. At the bottom right, there are "OK" and "Cancel" buttons.

SOP Class UID	1.2.840.10008.5.1.4.1.1.6	Key Image	0
SOP Instance UID	1.2.528.1.1001.1.19961112013.200.4.1996	Patient Orientation	NA
Number	1	Time	7:38:00 AM
Overlay Number	NA	Date	4/6/1993
Lut Number	NA	Series UID	1.2.528.1.1001.1.951018001.200.3.199601
Curve Number	NA		

OK Cancel

FIG. 15

128

Series Info

Modality	US	Device Serial Number	NA
Number		Software Version	NA
Description	NA	Study UID	1.2.528.1.1001.1.960113006.200.2.1
Operator Identification	NA	DICOM UID	1.2.528.1.1001.1.951018001.200.3.1
Performing Physician	NA	Institution Name	
Secondary Capture Device ID	NA	Study Status ID	
Scheduled Station Name	NA	Manufacturer	
Requested Contrast Bolus Agent	NA	Manufacturer Model Name	NA
Contrast Bolus Agent	NA		
Body Part	head		

OK Cancel

FIG. 16

130

Study Info		Patient Info at Time of Study	
DICOM UID:	1.2.528.1.1001.1.960113006.200.2.19960305.183737	Patient ID:	P0104
ID:	S023	Sex:	M
Time:	7:38:00 AM	Age:	NA
Date:	4/6/1993	Weight:	NA
Referring Physician:	Dr. Hu	Size:	NA
Reading Physician:	NA	Occupation:	NA
Accession Number:		Telephone Number:	NA
Description:	Brain	Address:	NA
Result UID:	NA	Country of Residence:	NA
Reason:	NA		
Modality In Study:	NA		
Admit Diagnosis:	NA		
Interpretation Author:	NA		

OK Cancel

FIG. 17

132

Form1

Event Type:

All Events

All Events

Image Database Mismatch

Invalid Tag Value Character Set (Image Level)

Invalid Tag Value Character Set (Patient Level)

Invalid Tag Value Character Set (Series Level)

Invalid Tag Value Character Set (Study Level)

Invalid Tag Value Length (Image Level)

Invalid Tag Value Length (Patient Level)

Invalid Tag Value Length (Series Level)

Invalid Tag Value Length (Study Level)

Patient Database Mismatch

Series Database Mismatch

Study Database Mismatch

☒ Group Like Events ☐ Include Reconciled

Apply Changes Cancel Changes

FIG. 18

134

Search Engine

Patient

General Info

ID: 136

Prefix:

First Name:

Middle Name:

Last Name:

Suffix:

Name Info

[-] Christensson^Jason

[-] S023 (Brain)

[-] US

1

2

3

4

5

6

7

8

9

10

11

12

13

14

Add New Study

Edit Existing Patent

Delete Patent

Merge Patent

View Patent

Search Advanced OK Cancel

FIG. 19

138

08002000	(0008,0020): Study Date. Decimal Tag: (0008,0032)
08000000	Tag Length: 8
3139393431323236	Tag Data: 19941226
08003000	(0008,0030): Study Time. Decimal Tag: (0008,0048)
0c000000	Tag Length: 12
3030303030302e3030303020	Tag Data: 000000.0000
08006000	(0008,0060): Modality. Decimal Tag: (0008,0096)
02000000	Tag Length: 2
4d52	Tag Data: MR
08007000	(0008,0070): Manufacturer. Decimal Tag: (0008,0112)
18000000	Tag Length: 24
5068696c697073204d65646963616c20	Tag Data: Philips Medical
53797374656d7320	Tag Data: Systems
08008000	(0008,0080): Institution Name. Decimal Tag:
(0008,0128)	
10000000	Tag Length: 16
444f4149204d454d4f5249414c204850	Tag Data: DOAI MEMORIAL HP
08003010	(0008,1030): Study Description. Decimal Tag:
(0008,4144)	
08000000	Tag Length: 8
432d5350494e4520	Tag Data: C-SPINE
10001000	(0010,0010): Patient's Name. Decimal Tag: (0016,0016)
12000000	Tag Length: 18
4d522f43205350494e452f5453452f32	Tag Data: MR/C SPINE/TSE/2
3536	Tag Data: 56
10003000	(0010,0030): Patient's Birth Date. Decimal Tag:
(0016,0048)	
08000000	Tag Length: 8
3139363430373331	Tag Data: 19640731
10004000	(0010,0040): Patient's Sex. Decimal Tag: (0016,0064)
02000000	Tag Length: 2
4d20	Tag Data: M

FIG. 20